



Entra Health Unveils MyGlucoHealth Clinical Point-of-Care System For Diabetes Specialists

Point-Of-Care Solution Provides A Clinic-based Testing Interface For Hospitals,
Clinics And Direct Inpatient Care

July 15, 2009 - San Diego, California – Entra Health Systems, an international applied healthcare technology company today unveiled the new MyGlucoHealth Clinical Point-of-Care System, a clinic-based diabetic testing solution for hospitals, clinics and other inpatient care environments. Clinical Point-of-Care is a hardware agnostic web-based interface that works with the MyGlucoHealth meter to allow blood glucose tests and data uploads using virtually any Bluetooth® or USB enabled tablet computer, PC or compatible PDA device.

“The Clinical Point-of-Care solution is designed to reduce information delays and improve accuracy through real-time, patient-direct documentation,” said Larry Mahar, Chief Technology Officer of Entra Health Systems. “MyGlucoHealth optimizes the clinical workflow by providing better, more accurate data to improve decision-making and collaboration, while decreasing the chance of transcription errors. This one step can reduce data documentation latency from hours to minutes. For clinical staff, it increases productivity and assists in better adherence to patient care protocols.”

Using the MyGlucoHealth meter in ‘point-of-care’ mode, a clinician can test multiple patients using a single meter. Readings can be separately uploaded and posted using built-in Bluetooth or USB to the MyGlucoHealth Physicians’ Portal or hospital EMR. On the Portal, readings for each patient are individually logged and color-coded based on patient criteria. Clinicians can log into the Portal independently to view a full patient dashboard, evaluate graphs, as well as set monitoring targets and alerts. The Portal will notify the clinician directly or alert the patient, family or caregivers if readings exceed thresholds. Permission to access Clinical Point-of-Care data can be granted to multiple clinicians to view a defined population of patient data.



Entra Health Systems

3111 Camino Del Rio North
Suite 101

San Diego, California 92108
US Office: 1(877)GLUCO-4-ME
info@entrahealthsystems.com
www.entrahealthsystems.com

In the United Kingdom
4A Blake Mews
Kew Gardens
Surrey TW9 3GA
United Kingdom
Phone: +44 (0) 113 815 5108